


Why Sanat Swing?

- ✓ High flex implant - 155° flexion
- ✓ Free combination of the implants
- ✓ CR, PS cemented and cementless, CC revision versions
- ✓ TiN coated version
- ✓ High intercondylar opening
- ✓ Innovative MIS instruments



Sanat Swing

TOTAL KNEE PROSTHESIS SYSTEM



- Cemented
- Cementless
- CC Revision

Innovative instruments for the quick and successful operation!




REV201703EN

Sanat Swing

TOTAL KNEE PROSTHESIS SYSTEM

COMBINATIONS



Sanat Swing primary knee prosthesis

- High flex (155° flexion) - anatomic movements ensured
- CR, PS versions available - wide range of applications
- Free combination of femoral and tibial implants
- Free combination of femoral and patellar implants
- TiN coating for metal allergic patients
- Cementless version

Instruments

- 4 trays CR +2 trays PS
- Step free varus/valgus setting
- Posterior referenced - perfect ligament tightness
- Extra- and intramedullary tibia targeting device
- Multifunction, minimally invasive instruments
- 35 minutes of operation time

CC - revision knee prosthesis

- Outstanding jump height in all movement phase
- Multiple offset options on both sides
- Tibia augments - huge variation possibilities
- Stems can be used on both sides
- Femoral augments for all cases - anterior augment is available
- Anatomic (following tibia shape) augmentation
- Augments can be fixed with screw
- Even increase in all femoral component dimension with size change
- Free combination of femoral and tibial components

Instruments

- 6 trays CC
- Same 3,5 mm screwdriver for all instruments
- Universal femoral sizing device
- Multi angle tibia resection is possible