


For precious hands



VORTEX Hand & Foot
Safe system for your success!


Contact

Address: 5 Faiskola str. Eger 3300, Hungary
E-mail: metal@sanatmetal.hu
Phone: +36 36 512 900



www.sanatmetal.com














Vortex Hand & Foot



Polyaxial angle stabilized system

www.sanatmetal.com

Vortex Hand & Foot Polyaxial angle stabilized system

Description	Field of Indication	Plate cross section (mm)	Image	Plate size	Plate length (mm)	Screw	
						Polyaxial	Cortical
1,5 XS straight plate	Hand	5x1,3		6H; 12H	38-77	1,5	1,5
1,5 XS T-plate	Hand	5x1,3		3-5H head 5-11H tail	38-77	1,5	1,5
1,5 XS L-plate	Hand	5x1,3		6H; 12H	38-77	1,5	1,5
1,5 XS parallel plate	Hand	5x1,3		4H; 6H; 8H	27-57	1,5	1,5
2,0 XS straight plate	Hand	6x1,8		6H; 12H	44-89	2,0	2,0
2,0 XS T-plate	Hand	6x1,8		3-5H head 5-11H tail	44-89	2,0	2,0
2,0 XS L-plate	Hand	6x1,8		6H; 12H	44-89	2,0	2,0
2,0 XS parallel plate	Hand	6x1,8		4H; 6H; 8H	33-69	2,0	2,0
Spider plate	Hand	17x1,8 20x1,8		small	-	2,7	2,7
				large	-		
Arthrodesis plate	Hand	8/12x3		4+5H	120	2,7	2,7
	Hand			4+4H	117,5	3,5	3,5
Mini straight plate	Hand, Foot	8x2,5		4H; 10H	47-101	2,7	2,7
Mini T-plate	Foot	7x2		3-5H head 5-11H tail	50-101	2,7	2,7
Mini L-plate	Foot	7x2		6H; 12H	50-101	2,7	2,7