



Sanatmetal 

SANAT R2



www.sanatmetal.com

SANAT SWING




Contact:

5 Faiskola str. Eger 3300, Hungary
+36 36 512-900
metal@sanatmetal.hu

www.sanatmetal.com

Sanatmetal 

PANNON



www.sanatmetal.com



HIP PROSTHESIS

ACETABULAR COMPONENTS

Cemented cup (UHMWPE/HXL-UHMWPE)	ConeTact cementless shell (Ti-TPS) and liner (UHMWPE/HXL-UHMWPE)		Bipolar head (Rex-steel/CoCr)
Cup N/L (5, 10 degrees) 42-62 mm Cup augment	Primer shell N/L (10, 20 degrees) 40-74 mm	REV shell N/L (10, 20 degrees) 40-74 mm Shell ring, flange	41-61 mm
			

FEMORAL HEADS (Rex-steel/CoCr/Ceramic)

22, 26, 28, 32, 36, 40 mm



STEMS

Cemented (Rex-steel)				Cementless (Ti)			
Sanat-E (Exeter) Centralizer, cement restrictor	Pannon (Müller)			Pannon (Müller)			Sanat R2
	C 1-10 sizes	C-LAT 1-10 sizes 125 degree	C-REV 3-10 sizes	CL-TPS 1-10 sizes	CL-HA 1-10 sizes	CL-LAT-HA 1-10 sizes 125 degree	CL-REV-TPS 3-10 sizes
0-5 sizes							
							

KNEE PROSTHESIS

FEMORAL COMPONENTS (CoCr)

Primer	Revision
Cemented, cementless CR/PS versions TIN coated A-H sizes	Cemented CCA version TIN coated A-H sizes Femoral anterior, distal, posterior augments Femoral extension stems and adaptors
	

TIBIAL INSERTS (UHMWPE)

Primer	Revision
CR/PS versions 1-8 sizes/10-22 mm	CC version 1-8 sizes/10-32 mm

TIBIAL COMPONENTS (CoCr, Ti)

Primer	Revision
Cemented, cementless TIN coated 1-8 sizes	Cemented TIN coated 1-8 sizes Tibial augments, extension stems and adaptors
	

PATELLA (UHMWPE)

Operated with cement, 26-42 mm